

Габаритные размеры

ТИП	DNA	DNM	a	g	f	x1	x2	n1	n2	n3	h1	h2	h3	m1	z1	s	b	i1	i	L	d	t	u	
40-65/2		65	40	180	268,5	380	439	409	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/3		65	40	250	268,5	380	509	479	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/4		65	40	320	268,5	380	579	549	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/5		65	40	390	268,5	380	649	619	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/6		65	40	460	268,5	380	719	689	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/7		65	40	530	268,5	380	789	759	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/8		65	40	600	268,5	380	859	829	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/9		65	40	670	268,5	380	929	899	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/10		65	40	740	268,5	380	999	969	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/11		65	40	810	268,5	380	1069	1039	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/12		65	40	880	268,5	380	1139	1109	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/13		65	40	950	268,5	380	1209	1179	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10
40-65/14		65	40	1020	268,5	380	1279	1249	280	212	160	180	180	180	60	60	15	60	120	109	80	32	35,3	10

